IKA HBR 4 digital
IKA HB 4 digital
IKA HB 4 basic

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60 010-1; EN 50 081; EN 50 082; EN 55 014 and EN 60 555.

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Guarantee

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Contents

	Page
Guarantee	2
Safety instructions	10
Intended use	10
Unpacking	10
Commissioning	10
Starting	11
Setting the safety circuit	11
Operating modes	11
Heating	12
Controlling the temperature of the medium	12
Stirring	12
Error messages	13
Maintenance and cleaning	13
Technical data	14
Associated standards and regulations	20
Accessories	20
List of spare parts HB 4 basic	22
Spare parts diagram HB 4 basic	23
List of spare parts HBR 4 digital / HB 4 digital	24
Spare parts diagram HBR 4 digital / HB 4 digital	25

Safety instructions

The safety temperature limit must be set at least 10K under the flashpoint of the tempering fluid used.

Never operate the heating baths without tempering fluid!

Filling and emptying when hot can lead to accidents, especially scalds.

Please observe the accident prevention regulations applicable for the job and place of work, wear protective goggles and gloves.

Take care when touching the external jacket, as it can heat up to over 80°C.

When working with temperature sensitive mediums please consider:

Upon continuous operation with high speed and a room temperature of 20°C a warming up of the housing up to 40°C is possible

When using an oil bath, only low-viscosity oils must be used, since using high-viscosity oils or solid grease can lead to overheating in localized areas, causing excessive pressure to build up in the bottom of the vessel.

When using an water bath, we recommend you use water with no mineral content.

Intended use

The IKA tempering baths HBR4 digital, HB4 digital and HB4 basic are laboratory devices and are suitable for directly tempering substances filled into the bath container.

Unpacking

Please unpack the equipment carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (post, rail or forwarder).

The delivery scope covers: an IKA tempering bath and operating instructions.

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

Starting

The device is switched on with the 2-pole ON/OFF mains switch (A). When switched on the lamp in the rocker switch lights up green.

Setting the safety circuit

IKA HBR4 digital and HB4 digital

When activating the HBR4 digital and HB4 digital devices an automatic self-test is carried out, all LEDs flash and all segments (**fig. 1**) are visible on the LCD display. After this, the LCD displays in **figure 2 to figure 5** are visible (standard factory setting).

When setting the safety circuit an upper temperature limit between 50°C and 210°C is set. The safety circuit can only be set in operating mode **"A"**.

- ◆ Switch on the device.
- ◆ Within 5 seconds of the display **fig. 6** appearing, keep the mode button **(E)** pressed down and set the safety temperature desired by simultaneously turning the control knob **(B)**.
- ◆ Release the mode button **(E)**. For two seconds the safety temperature set and safe OK **fig. 7** is displayed; the safety temperature set is saved.

The safety circuit is automatically checked each time the device is restarted.

IKA HB4 basic

A safety temperature limiter takes precedence over the main controller which works independently from the closed loop control circuit. The safety temperature is set with a small screw driver on the setting axle **(K)**. The end of the slot opposite the flattening is used as a setting mark.

If the heating bath is used as a water bath, the safety temperature limiter is to be set to the symbol **H₂Omax**.

When the safety temperature is reached the device switches off permanently. Only once the heating bath has cooled down somewhat can it be restarted by pressing in the red safety button **(L)** for example using a small screwdriver.

Operating modes

IKA HBR4 digital and HB4 digital

The device can be operated in three different operating modes (display see **fig. 8 to fig. 10**). To toggle from one operating mode to the other keep the mode button **(E)** depressed when switching on the device (min. 5 secs). Sequence A-B-C-A-B-C-A- etc.

Operating mode "A" This operating mode is the factory setting of the device. When switching on the device the heating and stirring functions are deactivated. The set values last set are saved and are used when the heating or stirring functions are switched on. The set values can be changed.

Operating mode "B" When switching on the device the state of the heating and stirring functions are taken over before the device was switched off last as well as the last values set. The set values can be changed.

Operating mode "C" When switching on the device the state of the heating and stirring functions are taken over before the device was switched off last as well as the last values set. The set values cannot be changed.

Heating

IKA HBR4 digital and HB4 digital

With the mode button **(E)** the correcting variable desired (set temperature) can be preselected. The heating bath temperature of the device is kept constant by the closed loop control circuit and is additionally monitored by the safety circuit. In the event of a disturbance in the closed loop control circuit, the heating bath is switched off permanently by the safety circuit. In the event of a fault in the closed loop control circuit or in the safety circuit, the LED "Temp." flashes green and yellow alternately. The LED next to the reset button lights up red.

The error is also displayed on the LCD (see chapter Error messages). The heating function cannot be started.

Controlling the temperature of the medium

IKA HB4 basic

The desired medium temperature can be set with the control knob **(B)** using the graduation on the front film. The yellow low voltage neon glow lamp **(J)** next to the control knob indicates the energy supply to the heating system of the device.

IKA HBR4 digital and HB4 digital

With the control knob **(B)** the medium temperature desired can be set between 0 and 200°C, not however higher than the safety temperature set. The value set can be read off on the digital display **(C)** (**fig. 11**). The heating function is switched on and off by pressing the sensor key "Temp. on/off" **(G)**. If this function is activated the green LED next to the lettering "Temp" lights up. The heating bath is heated up to the temperature set. On the LCD display the set temperature, actual temperature and the PV symbol are displayed. The set temperature and actual temperature relate to the medium.

During the heating phase the colour of the LED changes between yellow and green. If the yellow LED lights up this means that the heating system is being supplied with energy. If the heating function is switched off using the sensor key and the medium has reached a temperature of above 50°C the yellow LED flashes (20% on 80% off). On the LCD display °C and HOT are displayed alternately (**fig. 12 and 13**).

The temperature of the heating element is limited by the safety temperature set. The medium temperature is controlled by a fuzzy logic controller. Thanks to this, the medium temperature recorded by a Pt500 is heated up as fast as possible to the temperature set without any overshooting. The fuzzy logic controller automatically adapts to the various heat capacities of the different heat transfer fluids. This guarantees an optimal temperature management without temperature drift and waviness.

Stirring

IKA HBR4 digital

With the mode button **(E)** the correcting variable desired (set speed) can be preselected. With the control knob **(B)** the desired speed of the stirrer drive can be set between 0 and 800 rpm in steps of 50 rpm. The value set can be read off on the digital display **(C)** (**fig. 14**). The stirring function is switched on and off by pressing the sensor key "Motor on/off" **(F)**. If this function is activated the green LED lights up next to the lettering "Mot." The motor runs gently up to the speed set. On the LCD display the set speed, actual speed and the PV symbol are displayed (**fig. 15**).

If both functions (heating and stirring) are switched off, the LCD display always indicates the set temperature set. When the heating function is switched on it is given priority on the LCD display. If the stirring function is started with the sensor key or called up using the mode key **(E)**, the system switches over to speed display for five seconds. With the mode button **(E)** the LCD display can be changed from temperature to speed and vice versa at any time. After five seconds the system switches back to the predominant mode.

Error messages

The following error messages can be displayed (**fig. 17**) and have the following causes described:

No.	Error	Remedy
Er 3	Internal device temperature above 76°C	Allow the device to cool down. Mains switch OFF/ON
Er 4	Motor speed deviates ± 300 rpm from the speed. Set speed greater than 350 rpm	Reduce load. Mains switch OFF/ON
Er 4	Stirring function turned on Actual speed = 0 rpm Indication after approx 30 sec	Mains switch OFF/ON Increase the speed or turn the stirring function
Er 20	Safety relay does not close	Safety circuit is set too low, check mains voltage (mains voltage must be larger than the nominal voltage - 10%). Mains switch OFF/ON
Er 21	Safety relay does not open	RESET key/service
Er 22	Disturbance before the last time the device was switched OFF or safety temperature is lower than the temperature of the medium	RESET button
Er 23	Adjustable safety circuit faulty	Mains switch OFF/ON / service
Er 24	Safety temperature set exceeded	Set safety temperature higher. Mains switch OFF/ON
Er 26	System dry, no medium in the bath container	Allow the device to cool down, then refill with medium. Mains switch OFF/ON
Er 27	Calibration error	Mains switch OFF/ON. When switching on, a calibration is automatically carried out.
Er 28	Sensor rupture	Mains switch OFF/ON
Er 29	Short circuit in the safety sensor or controller and safety sensor mixed up	Service

Er 30 | Short circuit in the controller sensor | Service

Maintenance and cleaning

The safety circuit must be checked over at least once a year by the user. To do this, fill the temperature control bath with 1l of water and in operating mode „A“ set the safety temperature to 100°C. Then set the target temperature to 80°C and start the heating function by pressing the Temp ON/OFF key. After the target temperature has been reached, switch the moderating bath off and back on again with the power switch and set the safety temperature to 70°C. The temperature of the medium will then be 10°C above the safety temperature, and the safety circuit will be activated. Error message Er 22 or Er 24 must be displayed.

The aparate is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by IKA to clean IKA devices. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from IKA. Wear the proper protective gloves during cleaning of the devices. Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using another than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

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CAN/CSA C22.2 (1010-1)

EN 55 011

Associated standards and regulations

Associated EU guidelines

EMV guidelines: 89/336/EC

Machine guidelines: 73/023/EWG

Construction in accordance with the following safety standards

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List of spare parts

HB 4 basic


Item	Designation
9	Pedestal
10	Braided conductor with fuse
14	Braided conductor with fuse
18	Heated pot
20	Rear part
21	Exterior jacket
22	Spacer sleeve
23	Insulation mat
24	Metal plate
25	Device plug
26	Cover
27	Cover
29	Rocker switch
30	Base
31	Control knob
33	Capillary tube controller
35	Front plate
37	Retaining plate
38	Upper section
39	Lamp
2001	Operating unit compl.
2002	Capillary tube limiter

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
IKA® Works, Inc.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

North America

IKA® Works, Inc.

2635 North Chase Pkwy SE
Wilmington, NC 28405-7419
USA
TEL. +1 800 733-3037
TEL. +1 910 452-7059
FAX +1 910 452-7693
E-mail: usa@ika.net

IKA® Works, (Asia) Sdn Bhd

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Asia - Australia

IKA® Works (Asia) Sdn Bhd

Lot 2 Jalan Indah 1/2
Taman Industri Rawang Indah
48000 Rawang
Selangor, Malaysia
TEL. +60 3 6093-3322
FAX +60 3 6093-3940
E-mail: ika@tm.net.my

IKA® Japan Y.K.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

Japan

IKA® Japan Y.K.

293-1 Kobayashi-cho
Yamato Koriyama Shi
639-1026 Japan
TEL. +81 74358-4611
FAX +81 74358-4612
E-mail: japan@ika.de

IKA® Works do Brasil Ltda.

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

South America

IKA® Works do Brasil Ltda.

Estrada do Guerenguê, 491
Taquara Jacarepagua, RJ
Rio de Janeiro
22713-000 Brasil
TEL. +55 21 2435-9600
FAX +55 21 2435-9601
E-mail: brasil@ika.net

IKA® Works Guangzhou

LABORATORY TECHNOLOGY
ANALYZING TECHNOLOGY
PROCESSING EQUIPMENT

China

IKA® Works Guangzhou

173-175 Friendship Road
Guangzhou Economic & Technological
Development Zone
Guangzhou 510730, P.R.CHINA
TEL. +86 20 8222-6772
FAX +86 20 8222-6776
E-mail: sales@ikagz.com.cn